

BATHROOM AIDS | USER GUIDE

rainbow care

# FS697 STAINLESS STEEL MOBILE COMMUNE



*For illustration purposes only*

# STAINLESS STEEL MOBILE COMMODE

## PRODUCT OVERVIEW



*For illustration purposes only*



*For illustration purposes only*

### FS697 Stainless Steel

The FS697 is an impressive 3-in-1 Multi-Functional Commode, featuring a Stainless Steel Lightweight Frame, and detachable parts for easy and convenient use.

- Flip Up Arm Rest
- Detachable Bucket
- Detachable Foot Rest


## STAINLESS STEEL MOBILE COMMODE

### PRODUCT FEATURES



#### LEG REST

Pull the catch release to remove detachable leg rest



#### ARM REST

Press button to release lock and flip up arm rest



Unlock



Lock

#### LOCKABLE WHEELS

Push up to unlock and down to lock

*For illustration purposes only*

### CARE INSTRUCTIONS

- When using mobility commode, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Commode.
- When getting in and out of the commode, please do not stand on the footplates.
- The footplates should be folded up beforehand, or removed.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or damage the Commode.
- Do not exceed the maximum load.

---

PRODUCT SPECIFICATIONS	FS697
Overall Length	82 cm or 32.3"
Overall Width	54 cm or 21.3"
Overall Height	98 cm or 38.6"
Seat Depth	44 cm or 17.3"
Seat Width	44 cm or 17.3"
Seat Height	50 cm or 19.6"
Front Wheel Diameter	12.5 cm or 5"
Back Wheel Diameter	12.5 cm or 5"
Safe Working Load	100 kg
Overall Weight	11 kg

---

Check our our full range of mobile commodes. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

