THUASNE®

LigaStep Pedi Walker


fr	Botte de marche pédiatrique à montants latéraux rigides.....	3
en	Pediatric walking boots with a rigid side shell.....	3
de	Pädiatrische Unterschenkel-Fußorthese mit steifen Seitenstützen.....	3
nl	Walker voor kinderen met laterale, rigide staanders.....	4
it	Stivale da deambulazione pediatrico con scocca laterale rigida.....	4
es	Bota de caminar pediátrica con montantes laterales rígidos.....	5
pt	Bota de caminhada pediátrica com suportes laterais rígidos.....	5
da	Pædiatrisk støttestøvle med stive stivere i siden.....	5
fi	Kokonaan jäykkärunkoinen lasten kävelysaapas.....	6
sv	Pediatriskt fotstöd med styvt lateralt skal.....	6
el	Παιδιατρικός νάρθηκας ποδοκνημικής με άκαμπτες πλευρικές ράβδους.....	7
cs	Dětská fixační ortéza s pevnými bočními výztuhami.....	7
pl	But ortopedyczny pediatriczny ze sztywnymi wzmocnieniami bocznymi.....	7
lv	Stingrs sānu balsta bērnu pastaigu zābaks.....	8
lt	Eiti padedantis įtvaras vaikams su standžiais šoniniais statramsčiais.....	8
et	Jäikade külgedega pediaatriline tugisaabas.....	9
sl	Pediatrični škorenj s trdimi stranskima opornikoma.....	9
sk	Detská členková dlaha s pevnými postrannými dlahami.....	9
hu	Gyerek járótalapas bokarögzítőt oldalmerevítővel.....	10
bg	Педиатрична ортеза за ходене тип „ботуш“ с твърди странични опори.....	10
ro	Cizmă pediatrică cu montanți laterali rigizi.....	11
ru	Детский сапог для ходьбы с жесткими боковыми стойками.....	11
hr	Pedijatrijska čizma s čvrstim kalupom.....	11
zh	儿童助行靴，带坚挺护腿束带.....	12
ar	حذاء للمشي للأطفال ذو دعائم جانبية صلبة.....	12

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Istruçiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - إرشادات الاستعمال

S	21 - 25
M	25 - 31
L	31 - 34

© Thuasne - 6172701 (2023-09)
© Studio Caterin

UK CA CE  THUASNE
120, rue Marius Aujan
92300 Levallois-Perret
France

www.thuasne.com

Ⓐ

Ⓑ

Ⓒ

Ⓓ

Ⓔ

Ⓕ

	S	21 - 25
	M	25 - 31
	L	31 - 34

fr	Taille - Pointure européenne
en	Size - European shoe size
de	Größe - Europäische Schuhgröße
nl	Maat - Europese schoenmaat
it	Taglia - Numero di scarpa europeo
es	Talla - Número de pie europeo
pt	Tamanho - Número do sapato europeu
da	Størrelse - Europæisk størrelse
fi	Koko - Kengän koko (EU)
sv	Storlek - Europeisk storlek
el	Μέγεθος - Ευρωπαϊκό νούμερο υποδημάτων
cs	Velikost - Velikost obuvi (Evropa)
pl	Rozmiar - Europejska numeracja obuwia
lv	Izmērs - Eiropas izmērs
lt	Dydis - Europos avalynės dydis
et	Suurus - Suurus Euroopa süsteemis
sl	Velikost - Velikost Evropa
sk	Velkost - Európska veľkosť topánok
hu	Méret - Európai cipőméret
bg	Размер - Европейски размер на обувката
ro	Mărime - Mărime europeană
ru	Размер - Европейский размер
hr	Veličina - Europski broj cipela
zh	尺码-鞋码 (欧盟)
ar	المقاس - المقاس الأوروبي

fr	Stabilisation	Immobilisation
en	Stabilisation	Immobilisation
de	Stabilisierung	Ruhigstellung
nl	Stabilisatie	Immobilisatie
it	Stabilizzazione	Immobilizzazione
es	Estabilización	Inmovilización
pt	Estabilização	Imobilização
da	Stabilisering	Støtte
fi	Stabiloi	Immobilisoi
sv	Stabilisering	Immobilisering
el	Σταθεροποίηση	Ακινητοποίηση
cs	Stabilizace	Znehybnění
pl	Stabilizacja	Unieruchomienie
lv	Stabilizēšana	Imobilizācija
lt	Stabilizavimas	Imobilizavimas
et	Stabiliseerimine	Liikumatuks muutmine
sl	Stabiliziranje	Imobilizacija
sk	Stabilizácia	Imobilizácia
hu	Stabilizálás	Rögzítés
bg	Стабилизация	Обездвижване
ro	Stabilizare	Imobilizare
ru	Стабилизация	Иммобилизация
hr	Stabilizacija	Imobilizacija
zh	稳定	固定
ar	الاستقرار	التثبيت

fr

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en

PEDIATRIC WALKING BOOT WITH RIGID SIDE UPRIGHTS

Description/Destination

Walking boot with rigid side uprights to immobilize and stabilize the lower third of the leg, ankle and foot.

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Rigid components: polyamide - glass fiber - rubber - ethylene vinyl acetate.

Textile components: polyamide - polyurethane - polyester.

Properties/Mode of action

Immobilizes and stabilizes the lower third of the leg, ankle and foot. 90° rigid uprights on both sides.

Designed to ensure minimal dorsiflexion of the metatarsal bone. Curved, anti-slip sole that offers patients a natural, comfortable walking motion.

Product can be worn on either foot.

Indications

Ankle/foot post-operative/post-traumatic immobilization.

Severe ankle/foot sprain.

Stable fracture of the lower 1/3 of the leg.

Conservative treatment of ankle/foot fracture.

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not apply the product in direct contact with broken skin.

Do not use in the event of known allergy to any of the components.

Do not use in the event of unstable fractures or fractures of the upper part of the tibia or fibula.

Do not use for patients weighing > 68 kg / 150 lbs.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

The choice of the device version is at the discretion of the healthcare professional.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

The systematic use of a long sock is recommended when wearing the device.

Ensure correct position and avoid folds when placing foam liner in the device.

Check the condition of the affected limb daily (with particular attention for patients with sensory deficit).

Do not bring the product into direct contact with greasy substances (ointment, cream, etc.).

Do not bring the product into direct contact with a dressing.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

Do not wear the device in a medical imaging machine.

Do not wear the device when driving a vehicle.

Do not expose the product to extreme temperatures (> 60°C).

It is recommended to adequately tighten the device to achieve immobilisation without restricting blood circulation. In the event of pain, tingling or numbness, loosen the straps or deflate the air pads and ask a healthcare professional for advice.

It is recommended that the venous thromboembolic risk be assessed before any immobilisation is indicated. Follow the advice of a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

If the device comes into contact with water, dry the textile part and wipe the rigid part well with a dry cloth.

Prolonged immobilisation can sometimes cause muscle weakness. When resuming walking, be careful of risks of falling.

In case of history of venous or lymphatic disorders, consult a healthcare professional.

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.

Possible risk of venous thrombosis.

Possible risk of back pain due to gait pattern modification.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

1. Take the walking boot out of its plastic cover, unfasten the self-fastening straps and open the liner.

2. Do not remove the protections from the side uprights. They must remain in place while you put the boot on.

3. Sit down and carefully place your lower leg and foot in the foam liner, checking that your heel is correctly positioned at the bottom of the boot .

Close the foam liner, applying it firmly around the lower leg and foot .

4. Make sure that your leg is properly positioned between the two rigid side uprights.

5. Once you have checked that the side uprights are properly positioned, you can remove the protections from the side uprights .

The self-fastening surface on the inner part of the uprights attaches to the foam liner.

6. Fasten the walking boot again, tightening the self-fastening straps from bottom  to top .

7. Try walking with the boot to check that it fits properly and that it is working correctly .

Taking off the walking boot: Loosen the straps and open the liner.

Garment care

Product can be washed in accordance with the conditions shown on this instruction leaflet and on the label. If the device comes into contact with water, dry the textile part and wipe the rigid part well with a dry cloth.

Rigid components:

Wash the rigid part with a moist cloth.

Textile components:

The soft part can be fully removed for washing.

Replace in the original location before next use.

Hand wash the liner in cold water. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not dry clean. Do not tumble-dry. Do not iron. Squeeze out excess water. Dry flat. Dry away from any direct heat source (radiator, sun, etc.).

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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Thuasne UK Ltd
Thuasne SAS – SIREN/RCS Nanterre 542 091 186 – capital 1 950 000 euros