BATHROOM AIDS | USER GUIDE

rainbow care

# **SHOWER CHAIR**





# **SHOWER CHAIR**

#### FT7305 ALUMINIUM PRODUCT OVERVIEW



For illustration purposes only

#### FT7600 ALUMINIUM PRODUCT OVERVIEW





# **SHOWER CHAIR**

#### FT7305 ALUMINIUM





FT7600 ALUMINIUM

For illustration purposes only

#### Shower Chair

Features a lightweight Aluminium frame with adjustable legs to suit every users' height preference with back rest.

#### PRODUCT FEATURES



### HEIGHT ADJUSTMENT

Press the button and pull / push to the level that you want



PRODUCT SPECIFICATIONS	FT7305	FT7600
Overall Length	52 cm or 20.5"	45 cm or 17.7"
Overall Width	44 cm or 17.3"	52 cm or 20.5"
Overall Height	82 cm - 92 cm 32.3" - 36.2"	77 cm - 87 cm 30.3" - 34.3"
Seat Depth	41 cm or 16.1"	34 cm or 13.4"
Seat Width	44 cm or 17.3"	40 cm or 15.7"
Seat Height	45 cm - 55 cm 17.7" - 21.7"	43 cm - 53 cm 16.9" - 20.9"
Safe Working Load	100 kg	100 kg
Overall Weight	4 kg	3 kg

#### CARE INSTRUCTIONS

- When using shower chair, please pay attention to people and/or items around to avoid collision.
- Do not exceed the maximum load.

Check our our full range of shower chair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.



© 6223 0904 • 6931 7425 
© 8501 1878 ≅ sales@rainbowcare.com.sg ⊕ rainbowcare.com.sg
I @rainbowcare @ rainbowcaresg 
Blk 3007 Ubi Road 1 #02-410 Singapore 408701