

BATHROOM AIDS | USER GUIDE

rainbow care

SHOWER CHAIR



For illustration purposes only

SHOWER CHAIR

FT7306 ALUMINIUM PRODUCT OVERVIEW



For illustration purposes only

FT7210G ALUMINIUM PRODUCT OVERVIEW



For illustration purposes only

SHOWER CHAIR

FT7306 ALUMINIUM



For illustration purposes only

FT7210G ALUMINIUM



Shower Chair

Features a lightweight Aluminium frame with adjustable legs to suit every users' height preference, seat handles to grasp on to, and a 12 hole seat with back rest.

PRODUCT FEATURES



HEIGHT ADJUSTMENT

Press the button and pull / push to the level that you want

For illustration purposes only

PRODUCT SPECIFICATIONS	FT7306	FT7210G
Overall Length	48 cm or 18.9"	55 cm or 21.7"
Overall Width	49.5 cm or 19.5"	50 cm or 19.7"
Overall Height	65 cm - 73 cm 25.6" - 28.7"	71 cm - 81 cm 28" - 31.9"
Seat Depth	30 cm or 11.8"	30.5 cm or 12"
Seat Width (Without Handles)	40 cm or 15.7"	39 cm or 15.3"
Seat Width (With Handles)	50 cm or 19.7"	49.5 cm or 19.4"
Seat Height	37 cm - 45 cm 15.7" - 19.7"	41 cm - 51 cm 16.1" - 20.0"
Safe Working Load	100 kg	90 kg
Overall Weight	3 kg	4 kg

CARE INSTRUCTIONS

- When using shower chair, please pay attention to people and/or items around to avoid collision.
- Do not exceed the maximum load.

Check our our full range of shower chair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

