

BATHROOM AIDS | USER GUIDE

rainbow care

# SHOWER CHAIR



*For illustration purposes only*

## SHOWER CHAIR

### FT7500 ALUMINIUM PRODUCT OVERVIEW



### FT7500 ALUMINIUM PRODUCT OVERVIEW



#### Shower Chair

An adjustable shower chair for use in the bathroom. With 4 anti slip legs, it provides a sturdy non slip support even in very wet environments.

## SHOWER CHAIR

---

### PRODUCT FEATURES



### HEIGHT ADJUSTMENT

Press the button and pull /  
push to the level that you  
want

*For illustration purposes only*

PRODUCT SPECIFICATIONS	FT7500
Overall Length	34 cm or 13.4 ”
Overall Width	32 cm or 12.6 ”
Overall Height	33 cm - 43 cm 13.0 ” - 16.9 ”
Seat Depth	32 cm or 12.6”
Seat Width	32 cm or 12.6”
Seat Height	40 cm or 15.7”
Safe Working Load	80 kg
Overall Weight	2 kg

## CARE INSTRUCTIONS

- When using shower chair, please pay attention to people and/or items around to avoid collision.
- Do not exceed the maximum load.

Check our our full range of shower chair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

