BATHROOM AIDS | USER GUIDE

rainbow care

STATIONARY COMMODE





STATIONARY COMMODE

328L2 ALUMINIUM PRODUCT OVERVIEW



FT4213STEEL PRODUCT OVERVIEW





STATIONARY COMMODE

328L2 ALUMINIUM



FT4213 STEEL



For illustration purposes only

Stationary Commode

Made of lightweight aluminium frame and functions as a 2 in 1 commode where it is used as a commode chair and a shower chair. Folable for storage.

PRODUCT FEATURES



HEIGHT ADJUSTMENT

Press the button and pull / push to the level that you want

REMOVE BUCKET

328L2 - Lift up the seat cover to remove FT4213 - Slide out to remove



PRODUCT SPECIFICATIONS	328L2	FT4213
Overall Length	57 cm or 22"	52 cm or 20.5"
Overall Width	56.5 cm or 22.4"	52.5 cm or 20.6"
Overall Height	87 cm - 91.5 cm 34.3" - 36"	68.5 cm - 78.5 cm 26.9" - 30.9"
Seat Depth	42 cm or 16.5"	42 cm or 16.5"
Seat Width	42 cm or 16.5"	42 cm or 16.5"
Seat Height	45 cm - 49.5 cm 17.7" - 19.5"	34.5 cm - 44.5 cm 13.6" - 17.5"
Safe Working Load	100 kg	80 kg
Overall Weight	5.5 kg	4.5 kg

CARE INSTRUCTIONS

- When using stationary commode, please pay attention to people and/or items around to avoid collision.
- Make sure the level of the height of each leg is the same before sitting down.
- Do not exceed the maximum load.

Check our our full range of stationary commodes. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.



© 6223 0904 • 6931 7425
© 8501 1878
∞ sales@rainbowcare.com.sg
⊕ rainbowcare.com.sg
☑ rainbowcaresg
② Blk 3007 Ubi Road 1 #02-410 Singapore 408701