BATHROOM AIDS | USER GUIDE rainbow Care STATIONARY COMMODE





STATIONARY COMMODE

329L2 ALUMINIUM PRODUCT OVERVIEW



FT4213STEEL PRODUCT OVERVIEW





STATIONARY COMMODE

329L2 ALUMINIUM



FT4213 STEEL



For illustration purposes only Stationary Commode

Made of lightweight aluminium frame and functions as a 2 in 1 commode where it is used as a commode chair and a shower chair. Folable for storage.

PRODUCT FEATURES



HEIGHT ADJUSTMENT

Press the button and pull / push to the level that you want

REMOVE BUCKET

329L2 - Lift up the seat cover to remove FT4213 - Slide out to remove



PRODUCT SPECIFICATIONS	329L2	FT4213
Overall Length	54 cm or 21.3″	52 cm or 20.5″
Overall Width	59 cm or 23.2″	52.5 cm or 20.6″
Overall Height	80 cm - 84.5 cm 31.5″ - 33.3″	68.5 cm - 78.5 cm 26.9″ - 30.9″
Seat Depth	40 cm or 15.7″	42 cm or 16.5″
Seat Width	42.5 cm or 16.7″	42 cm or 16.5″
Seat Height	43 cm - 47.5 cm 16.9″ - 18.7″	34.5 cm - 44.5 cm 13.6″ - 17.5″
Safe Working Load	130 kg	80 kg
Overall Weight	5.5 kg	4.5 kg

CARE INSTRUCTIONS

- When using stationary commode, please pay attention to people and/or items around to avoid collision.
- Make sure the level of the height of each leg is the same before sitting down.
- Do not exceed the maximum load.

Check our our full range of stationary commodes. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.



© 6223 0904 • 6931 7425
 © 8501 1878
 ∞ sales@rainbowcare.com.sg
 ⊕ rainbowcare.com.sg
 ☑ rainbowcaresg
 ② Blk 3007 Ubi Road 1 #02-410 Singapore 408701