

ZINGER ELECTRICAL WHEELCHAIR



ZINGER ELECTRICAL WHEELCHAIR

PRODUCT OVERVIEW



Zinger Electrical Wheelchair



Zinger sets new standards with its ultra simple, reliable traction drive and intuitive steering system, giving it amazing manoeuvrability. Its unique design creates a revolutionary instant folding frame.

ELECTRICAL WHEELCHAIR

PRODUCT FEATURES




INSERT BATTERY

Always insert one side first by aligning the catch



TO REMOVE THE BATTERY

Turn the key clockwise to unlock and pull up the battery



For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

ZINGER ELECTRICAL WHEELCHAIR

PRODUCT SPECIFICATIONS

Overall Length	79 cm or 31"
Overall Width (Wheel to Wheel)	64.5 cm or 25"
Overall Height	93.5 cm or 37"
Folded Length	95.5 cm or 38"
Folded Width	64.5 cm or 25"
Folded Height	25 cm or 10"
Seat Depth	38.1 cm or 15"
Seat Width	42 cm or 16.5"
Seat Height	48.3 cm or 19"
Arm Rest Height (Floor to Arm Rest)	63.5 cm or 25"
Back Rest Height	40.5 cm or 16"
Front Wheel Diameter	15 cm or 6"
Back Wheel Diameter	22.8 cm or 9"
Safe Working Load	113 kg
Overall Weight	19 kg

Check our our full range of electrical wheelchair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

